

Infrared Thermometer, Mintti Breeze



Mintti Breeze infrared forehead and ear thermometer integrates with the latest built-in smart chip and AI algorithm. It can measure temperature with the highest degree of medical accuracy for all ages. Whether you are at home or work as a physician, Mintti Breeze allows you to monitor your temperature easily.

To meet the users' different needs and circumstances, Mintti Breeze can measure temperatures with contact or contact-free methods. The contact-free method with the silent mode is ideal for infants without waking them up while taking the temperature. You can also easily switch to ear reading to ensure more accurate reading on older children and adults.

With the world's first front-facing LED display, you don't have to turn Mintti Breeze over to check the reading as the display is right in front of you. It can make taking temperature so convenient. The need for self-measuring becomes more comfortable and much easier.

SPECIFICATION:

Product Name	Mintti Breeze Infrared Thermometer
Product Size	161.5mm *35.8mm *55.9mm(length x width x height)
Product Weight	62 g (without battery)
Measurement Temperature	Ear Mode 34.0 °C~42.4 °C (93.2 °F~108.3 °F)
	Forehead Mode 34.0 °C~42.4 °C (93.2 °F~108.3 °F)
Measurement Accuracy	Ear Mode ± 0.2 °C/0.4 °F
	Forehead Mode ± 0.2 °C/0.4 °F
Bluetooth	4.0 Enabled, Smart APP
Memory Storage	Up to 32 Readings
Storage Tem. & Hum.	-20 ~ 55 °C(-4 °F to 131 °F) & $\leq 90\%$ moisture noncondensing
Power Supply	DC1.5V X 2AAA
Power Consumption	When off ≤ 10 UW When measuring ≤ 200 MW
Auto Power-off	30 seconds
Product Accessories	AAA batteries x2 and User Manual x1 Maintenance Service Card x 1 and QC Certificate x 1

FEATURES:

Fast & Accurate Measuring

Intelligent Performance & User Friendly

Combining the cutting-edge sensor with a smart chip and optimized AI algorithms, Mintti Breeze can prevent false reading and offer the highest degree of accuracy within one second.

Easy-to-Read Front Facing LED Display

With the world's first front-facing display, it makes self-measuring so easy. There is no need to turn it over to check the temperature as it is right in front of you.

Preventing Infection with the Contact-free Forehead Measuring Mode

With the contact-free forehead measuring mode, it can prevent infection and cross-contamination and minimize the disease from spreading among users.

It features with a built-in fever alarm to prevent guesswork and an auto shutdown for power saving. With one-touch operation, it is simple to use even for children and elderly people.

More Durable with Long-lasting Battery

The long-lasting battery with approximately 3000 readings per year and 200 minutes for constant use makes it more reliable and durable to use.

Ergonomic and Elegant Patented Design

The ergonomic design of Mintti Breeze makes taking temperature more comfortable and much easier. With a unique upright patented design, it looks more like an artwork rather than a thermometer.