

Infrared Thermometer, Mintti Breeze



The Mintti Breeze infrared forehead and ear thermometer is equipped with the latest integrated smart chip and AI algorithm, ensuring the highest degree of medical accuracy for temperature measurement across all age groups. Whether you're at home or working as a physician, Mintti Breeze allows you to easily monitor your temperature.

To cater to various user needs and situations, Mintti Breeze offers both contact and contact-free temperature measurement methods. The contact-free method, featuring a silent mode, is ideal for measuring infants' temperatures without disturbing their sleep. Additionally, you can effortlessly switch to ear reading for more precise temperature readings in older children and adults.

Featuring the world's first front-facing LED display, Mintti Breeze eliminates the need to flip it over to check the reading, as the display is conveniently positioned in front of you. This

design greatly enhances the convenience of temperature measurement, making self-monitoring more comfortable and effortless.

SPECIFICATION:

Product Name	Mintti Breeze Infrared Thermometer
Product Size	161.5mm *35.8mm *55.9mm(length x width x height)
Product Weight	62 g (without battery)
Measurement Temperature	Ear Mode 34.0°C ~42.4°C (93.2°F ~ 108.3°F)
	Forehead Mode 34.0°C ~42.4°C (93.2°F ~ 108.3°F)
Measurement Accuracy	Ear Mode ±0.2°C/0.4°F
	Forehead Mode ±0.2°C/0.4°F
Bluetooth	4.0 Enabled, Smart APP
Memory Storage	Up to 32 Readings
Storage Tem. & Hum.	-20°~ 55°C(-4°F to 131°F) & ≤90% moisture noncondensing
Power Supply	DC1.5V X 2AAA
Power Consumption	When off ≤ 10 UW When measuring ≤ 200 MW
Auto Power-off	30 seconds
Product Accessories	AAA batteries x2 and User Manual x1 Maintenance Service Card x 1 and QC Certificate x 1

FEATURES:

Fast & Accurate Measuring

Combining cutting-edge sensors with a smart chip and optimized AI algorithms, Mintti Breeze can prevent false readings and offer the highest degree of accuracy within one second.

Easy-to-Read Front Facing LED Display

With the world's first front-facing display, it makes self-measuring so easy. There is no need to turn it over to check the temperature as it is right in front of you.

Preventing Infection with the Contact-free Forehead Measuring Mode

The contact-free forehead measuring mode can prevent infection and cross-contamination, minimizing the spread of disease among users.

Intelligent Performance & User Friendly

It features a built-in fever alarm to prevent guesswork and an auto shutdown for power saving. With one-touch operation, it is simple to use even for children and elderly people.

More Durable with Long-lasting Battery

The long-lasting battery with approximately 3000 readings per year and 200 minutes of constant use makes it more reliable and durable to use.

Ergonomic and Elegant Patented Design

The ergonomic design of Mintti Breeze makes taking temperature more comfortable and much easier. Its unique upright patented design makes it look more like artwork rather than a thermometer.